



JPD-100B+

Fetal Doppler

Features

- ▲ Safe: The ultrasound is lower than the national ultrasonic power standard by over 50%.
- ▲ Portable: Compact and portable design provides a convenient use and carrying.
- ▲ LCD display: Large backlight LCD displays the fetal heart rate.
- ▲ Loudspeaker: Built-in speaker gives clear and loud voice.
- ▲ Interchangeable waterproof probe
- ▲ Automatic shutdown function in case of no work for one minute.
- ▲ Fully charged battery enables to work continuously for over 10 hours.
- ▲ High sensibility can detect the pulse of 9-week fetus.
- ▲ Three calculation modes to calculate fetal heart rate: instant, 10 second average and manual

Specification

- ▲ Display: LCD display
- ▲ Heart Rate Range: 50-210BPM
- ▲ Resolution: ± 1 BPM
- ▲ Audio output power: 0.5w
- ▲ Auto Power-off: 1 minute with no signal or operation
- ▲ Size: 135mm*84mm*33mm
- ▲ Weight: Main body 0.43kg



Shenzhen Jumper Medical Equipment Co., Ltd.

Add: 5th Floor, Building No. 34, Baoyuan Industrial Zone, Xixiang Street, Baoan District, Shenzhen

Tel: 0086-755-26696279

FAX: 0086-755-26852025

Email: info@jumper-medical.com

Web: www.jumper-medical.co