

JPD-300P

FETAL MONITOR



Jumper[®]

JPD-300P

FETAL MONITOR

[Product introduction]

- ▲ The monitor of lightweight and space-saving.
- ▲ It can meet the need of whole stages of labor(including antepartum and intrapartum monitoring for fetus and mother).
- ▲ It also can be connected to the obstetrical central station



[Options]

- ▲ Maternal Spo2,HR,NIBP and temperature measurement.
- ▲ Maternal ECG and PR measurement.
- ▲ Twins monitor interface,Twins FHR measurement.
- ▲ Build-in network/PC/connecting interface
- ▲ Fas(fetal acoustic stimulator)



[Features]

- ▲ Light weight,compact design,space-saving,easy operation
- ▲ Interface of parameters of FHR,TOCO,FM,with digital and curve display.
- ▲ The setup of the system can be easy done and automatically stored.
- ▲ Alarm system.
- ▲ 12-inch TFT color Folding up to 90 degree.
- ▲ Build-in152 mm long life thermal printer.
- ▲ DSP technology,real time identify and measurement for FHR,accurate and reliable results.
- ▲ Special high sensitive water proof probe
- ▲ More than 24 hours data storage which can be playback and printed
- ▲ Special patient event and clinical event marking button.
- ▲ Low ultrasound power,to fetus safety
- ▲ Wide range of applying voltage(100-240V)
- ▲ Automatic Fetal Movement Detection.



SHENZHEN JUMPER MEDICAL EQUIPMENT CO.,LTD.

5th Floor, Building No.34, Baoyuan Industrial Zone, Xixiang Street, Baoan District, Shenzhen Zip Code: 518102
Tel: 86-755-26696279 Website : www.jumper-medical.com email : info@jumper-medical.com