

JPD-800A

Patient Monitor [7 inch]



[Feature]

7 inch TFT display

Unique Human voice alarm

Maximum 720 hour graphic and tabular trends

Maximum 12 hours of all waveforms

Support USB data output

System software can be upgraded through network

Wire and wireless network





[Specifications]

ECG

lead mode;I、II、III、AVR、AVL、AVF、V
Gain:X0.25、X0.5、X1、X2
Heart Rate:15-300BPM(Adult) 15-350BPM(Neonatal)
Resolution:1 BPM
Accuracy:+-1%
ST measurement Range:-2.0-+2.0mV
Accuracy:+-10%
Pacemaker:Yes
Sweep Speed:12.5 mm/s,25mm/s,50mm/s
Bandwidth:Diagnostic:0.05--130Hz
Monitor: 0.5--40Hz
Surgery: 1-20Hz

SPO2

Measure Range:0-100%
Resolution:1%
Accuracy:2%(70-100%)
Pulse Rate:20-300BPM
Resolution:1BPM
Accuracy:+-3BPM

RESP:

Measure Range:0-150BrPM
Resolution:1BrPM
Accuracy:+-2%

Temp:

Measure Range:0-50C
Resolution:0.1C
Accuracy:+-0.1C

NIBP:

Work mode:Manual,Automatic,Continuous
Unit:mmHg,Kpa
Measure Range:
Adult: SYS: 40-280mmHg
DIA: 10-220mmHg
MEAN:20-240mmHg

Pediatric:

SYS:40-220mmHg
DIA: 10-160mmHg
MEAN:20-170mmHg

Neonatal:

SYS:40-135mmHg
DIA:10-100mmHg
MEAN:20-110mmHg

Resolution:1mmHg
Accuracy:+-5mmHg



SHENZHEN JUMPER MEDICAL EQUIPMENT CO.,LTD.

Add: 5th Floor, Building No.34,Baoyuan Industrial Zone,Xixiang Street, Baoan District,Shenzhen
Tel: +86-755-26696279 Fax: +86-755-26852025 Email:info@jumper-medical.com
Web: www.jumper-medical.com